


Anti-Xpd/Ercc2 Polyclonal Antibody For Wb, If/Icc, Elisa - P18074

Item Number: CM0024817



Introduction

Rabbit polyclonal antibody against human XPD/ERCC2 (TFIIH helicase). Validated for WB, IF/ICC, and ELISA applications. Cross-reacts with human and rat. Ideal for DNA damage repair and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	XPD/ERCC2 Rabbit pAb
Remarks/Alias	EM9; TTD; XPD; TTD1; COFS2; TFIIH; XPD/ERCC2
Species	Human
Gene ID (Human)	2068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human XPD/ERCC2 (NP_001124339.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P18074
Protein Weight	87kDa
Shipping	Ice bag