

Anti-Xpnpep2 Rabbit Polyclonal Antibody For Wb And Elisa - O43895

Item Number: CM0013819

Introduction

Rabbit polyclonal antibody targeting XPNPEP2 (X-Pro aminopeptidase 2). Suitable for WB and ELISA, cross-reacts with human, mouse, and rat. Generated against recombinant ectodomain (aa 22-240). Useful for bradykinin metabolism research.

[Learn More](#)

Parameter	Value
Product Name	XPNPEP2 Rabbit pAb
Remarks/Alias	APP2; AEACEI; XPNPEP2
Species	Human
GeneID (Human)	7512
GeneID	7512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-240 of human XPNPEP2 (NP_003390.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43895
Protein Weight	76kDa
Shipping	Ice bag