

Anti-Xpnpep3 Polyclonal Antibody (Rabbit) For Wb, Elisa - Q9Nqh7

Item Number: CM0020429



Introduction

Rabbit polyclonal antibody against human XPNPEP3 (Xaa-Pro aminopeptidase 3), validated for WB and ELISA. Detects human and mouse XPNPEP3; suitable for mitochondrial peptidase and TNF signaling research.

[Learn More](#)

Parameter	Value
Product Name	XPNPEP3 Rabbit pAb
Remarks/Alias	XPNPEP3; APP3; ICP55; NPHPL1; X-prolyl aminopeptidase 3
Species	Human
GeneID (Human)	63929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-230 of human XPNPEP3 (NP_071381.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQH7
Protein Weight	31kDa/32kDa/48kDa/54kDa/57kDa
Shipping	Ice bag