

Anti-Xpnpep3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nqh7

Item Number: CM0014980



Introduction

Rabbit polyclonal antibody against XPNPEP3 (APP3). Validated for WB, IF/ICC, ELISA; detects human, mouse, rat. For research on mitochondrial aminopeptidase, TNF signaling, and nephronophthisis.

[Learn More](#)

Parameter	Value
Product Name	XPNPEP3 Rabbit pAb
Remarks/Alias	APP3; ICP55; NPHPL1; XPNPEP3
Species	Human
Gene ID (Human)	63929
Gene ID	63929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human XPNPEP3 (NP_071381.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQH7
Protein Weight	57kDa
Shipping	Ice bag