

Anti-Xpo6 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Qu8

Item Number: CM0024204

Introduction

Rabbit polyclonal antibody against human XPO6 (Exportin-6 / EXP6 / RANBP20), validated for Western blot (WB) and ELISA. Detects human and mouse XPO6. Suitable for nucleocytoplasmic transport and actin dynamics research.

[Learn More](#)

Parameter	Value
Product Name	XPO6 Rabbit pAb
Remarks/Alias	XPO6; EXP6; RANBP20; exportin-6
Species	Human
Gene ID (Human)	5159
Gene ID	23214
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 116-216 of human XPO6 (NP_055986.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96QU8
Protein Weight	127kDa/128kDa
Shipping	Ice bag