

Anti-Xpot Rabbit Polyclonal Antibody For Wb, Elisa - O43592

Item Number: CM0014396

Introduction

Rabbit polyclonal antibody targeting exportin-T (XPOT/XPO3) for Western blot and ELISA. Reacts with human, mouse, and rat XPOT. Ideal for studying tRNA nuclear export transport mechanisms.

[Learn More](#)

Parameter	Value
Product Name	XPOT Rabbit pAb
Remarks/Alias	XPO3; XPOT
Species	Human
GeneID (Human)	11260
GeneID	11260
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human XPOT (NP_009166.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43592
Protein Weight	110kDa
Shipping	Ice bag