

Anti-Xrcc1 Rabbit Monoclonal Antibody For Wb, Elisa - P18887

Item Number: CM0011386



Introduction

High-quality rabbit monoclonal antibody against XRCC1, validated for WB and ELISA. Detects human, mouse, and rat XRCC1, a key scaffold protein in DNA single-strand break repair.

[Learn More](#)

Parameter	Value
Product Name	XRCC1 Rabbit mAb
Remarks/Alias	RCC; SCAR26; XRCC1
Species	Human
Gene ID (Human)	7515
Gene ID	7515
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18887
Protein Weight	69kDa
Shipping	Ice bag