

Anti-Xrcc2 Polyclonal Antibody For Wb, Ihc-P, Elisa - O43543

Item Number: CM0022321



Introduction

Rabbit polyclonal antibody targeting human XRCC2, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for DNA repair and homologous recombination studies.

[Learn More](#)

Parameter	Value
Product Name	XRCC2 Rabbit pAb
Remarks/Alias	FANCU; POF17; SPGF50; XRCC2
Species	Human
GeneID (Human)	7516
GeneID	7516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-83 of human XRCC2 (NP_005422.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43543
Protein Weight	32 kDa
Shipping	Ice bag