


Anti-Y14/Rbm8A Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y5S9

Item Number: CM0002514



Introduction

Rabbit monoclonal antibody against human Y14/RBM8A, validated for WB and ELISA. Reacts with human, mouse, rat. Ideal for studying exon junction complex (EJC), pre-mRNA splicing, and nonsense-mediated mRNA decay pathways.

[Learn More](#)

Parameter	Value
Product Name	Y14/RBM8A Rabbit mAb
Remarks/Alias	TAR; Y14; RBM8; ZNRP; RBM8B; ZRNP1; BOV-1A; BOV-1B; BOV-1C; MDS014; DEL1q21.1; C1DELq21.1; Y14/RBM8A
Species	Human
GeneID (Human)	9939
GeneID	9939
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Y14/RBM8A (Q9Y5S9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5S9
Protein Weight	20kDa
Shipping	Ice bag