

Anti-Yap1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P46937

Item Number: CM0003472



Introduction

Anti-YAP1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Key transcription coactivator in Hippo pathway, organ size control, and tumor suppression. Recombinant protein immunogen.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] YAP1 Rabbit mAb
Remarks/Alias	YAP, YKI, COB1, YAP2, YAP-1, YAP65, P1
Species	Human
GeneID (Human)	10413
GeneID	10413
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-504 of human YAP1 (NP_001123617.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46937
Protein Weight	36 kDa-54 kDa
Shipping	Ice bag