


## Anti-Yap1 Rabbit Monoclonal Antibody For Wb - P46937

Item Number: CM0006312



### Introduction

YAP1 rabbit monoclonal antibody (KO validated) for WB, IHC-P, IF/ICC, IP, ELISA, ChIP. Detects human, mouse, rat YAP1. Targets amino acids 155-504. Ideal for Hippo signaling, organ size control, and tumor suppression studies.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] YAP1 Rabbit mAb                                                                                       |
| Remarks/Alias    | YAP; YKI; COB1; YAP2; YAP-1; YAP65; P1                                                                               |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 10413                                                                                                                |
| GeneID           | 10413                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 155-504 of human YAP1 (NP_001123617.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Monoclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA, ChIP                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | P46937                                                                                                               |
| Protein Weight   | 54kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |