


Anti-Yap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P46937

Item Number: CM0014490



Introduction

Rabbit polyclonal antibody targeting human YAP1, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Recognizes the Hippo pathway transcriptional coactivator YAP1 (P46937). Ideal for cell signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	YAP1 Rabbit pAb
Remarks/Alias	YAP; YKI; COB1; YAP2; YAP-1; YAP65; YAP1
Species	Human
GeneID (Human)	10413
GeneID	10413
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-504 of human YAP1 (NP_001123617.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46937
Protein Weight	54kDa
Shipping	Ice bag