


## Anti-Yap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P46937

Item Number: CM0018442



### Introduction

Rabbit polyclonal antibody to human YAP1, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Ideal for Hippo signaling and cancer research.

[Learn More](#)

| Specification    | Details                                                                                                              |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | YAP1 Rabbit pAb                                                                                                      |
| Remarks/Alias    | YAP; YKI; COB1; YAP2; YAP-1; YAP65; YAP1                                                                             |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 10413                                                                                                                |
| Gene ID          | 10413                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 155-504 of human YAP1 (NP_001123617.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                         |
| SWISS            | P46937                                                                                                               |
| Protein Weight   | 54kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |