


Anti-Yars Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P54577

Item Number: CM0028740



Introduction

Tyrosine--tRNA ligase (YARS) Rabbit Polyclonal Antibody validated for WB, IP, ELISA in human, mouse, and rat. Suitable for research on aminoacyl-tRNA synthetases, protein synthesis, and PARP1 regulation.

[Learn More](#)

Parameter	Value
Product Name	YARS Rabbit pAb
Remarks/Alias	YRS; YTS; YARS; TYRRS; CMTDIC; IMNEPD2
Species	Human
GeneID (Human)	8565
GeneID	8565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human YARS (NP_003671.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54577
Protein Weight	59kDa
Shipping	Ice bag