


Anti-Yb-1/Ybx1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P67809

Item Number: CM0011873



Introduction

High-quality rabbit monoclonal antibody against YB-1/YBX1. Suitable for WB, IP, IF, IHC, ELISA. Detects human, mouse, rat YB-1. Ideal for RNA processing, transcriptional regulation, and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	YB-1/YBX1 Rabbit mAb
Remarks/Alias	YB1; BP-8; CSDB; DBPB; YB-1; CBF-A; CSDA2; EFI-A; NSEP1; NSEP-1; MDR-NF1; YB-1/YBX1
Species	Human
Gene ID (Human)	4904
Gene ID	4904
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 225-324 of human YB-1/YBX1 (P67809).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67809
Protein Weight	36 kDa
Shipping	Ice bag