


Anti-Yb-1/Ybx1 Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P67809

Item Number: CM0025845



Introduction

YB-1/YBX1 Rabbit Polyclonal Antibody validated for WB, IF-P, IHC-P, ELISA; reacts with human, mouse, rat YB-1/YBX1 protein (36kDa). Ideal for studying mRNA stability, splicing, and exosome biology.

[Learn More](#)

Parameter	Value
Product Name	YB-1/YBX1 Rabbit pAb
Remarks/Alias	YB1; BP-8; CSDB; DBPB; YB-1; CBF-A; CSDA2; EFI-A; NSEP1; NSEP-1; MDR-NF1; YB-1/YBX1
Species	Human
GeneID (Human)	4904
GeneID	4904
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human YB-1/YBX1 (NP_004550.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67809
Protein Weight	36kDa
Shipping	Ice bag