


Anti-Yb-1/Ybx1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P67809

Item Number: CM0030076



Introduction

High-quality rabbit polyclonal antibody targeting human YB-1/YBX1, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for transcription, RNA biology, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	YB-1/YBX1 Rabbit pAb
Remarks/Alias	YB1; BP-8; CSDB; DBPB; YB-1; CBF-A; CSDA2; EFI-A; NSEP1; NSEP-1; MDR-NF1; YB-1/YBX1
Species	Human
GeneID (Human)	4904
GeneID	4904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-55 of human YB-1/YBX1 (NP_004550.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67809
Protein Weight	36kDa
Shipping	Ice bag