

Anti-Ybey Rabbit Polyclonal Antibody For Wb And Elisa - P58557

Item Number: CM0022768

Introduction

High-quality rabbit polyclonal antibody against human endoribonuclease YBEY, validated for Western blot and ELISA with cross-reactivity to rat, ideal for rRNA maturation and nucleolar research.

[Learn More](#)

Parameter	Value
Product Name	YBEY Rabbit pAb
Remarks/Alias	C21orf57
Species	Human
GeneID	54059
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 68-167 of human YBEY (NP_478061.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P58557
Protein Weight	19kDa
Shipping	Ice bag