


Anti-Yeats4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95619

Item Number: CM0026623



Introduction

Rabbit polyclonal antibody against human YEATS4 (GAS41), a chromatin reader in the NuA4 complex. Suitable for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat YEATS4.

[Learn More](#)

Parameter	Value
Product Name	YEATS4 Rabbit pAb
Remarks/Alias	YAF9; GAS41; NUBI-1; 4930573H17Rik; B230215M10Rik; YEATS4
Species	Human
GeneID (human)	8089
GeneID	8089
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human YEATS4 (NP_006521.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95619
Protein Weight	26kDa
Shipping	Ice bag