

Anti-Yipf6 Monoclonal Antibody For Wb, Elisa - Q8Br70

Item Number: CM0011723

Introduction

Anti-Yipf6 rabbit monoclonal antibody validated for WB and ELISA against mouse YIPF6 protein, a Golgi apparatus membrane protein involved in YIPF1 and YIPF2 stability. Suitable for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Yipf6 Rabbit mAb
Remarks/Alias	A430107J06Rik
Species	Mouse
GeneID (Human)	286451
GeneID	77929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-83 of mouse Yipf6 (NP_997516.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8BR70
Protein Weight	26kDa
Shipping	Ice bag