

Anti-Yipf6 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Br70

Item Number: CM0024683

Introduction

Rabbit polyclonal antibody targeting mouse YIPF6 (YIP1 family member 6). Validated for WB and ELISA with cross-reactivity to mouse and rat. Explores Golgi membrane protein function and YIP1 family interactions.

[Learn More](#)

Parameter	Value
Product Name	Yipf6 Rabbit pAb
Remarks/Alias	A430107J06Rik; Yipf6
Species	Mouse
GeneID (Human)	286451
GeneID	77929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-83 of mouse Yipf6 (NP_997516.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8BR70
Protein Weight	26kDa
Shipping	Ice bag