

Anti-Yipf7 Polyclonal Antibody For Wb And Elisa - Q8N8F6

Item Number: CM0020363

Introduction

Rabbit polyclonal antibody targeting human YIPF7 (FinGER9, YIP1 family member 7). Suitable for Western blot and ELISA, with cross-reactivity to rat. Ideal for Golgi/ER trafficking research.

[Learn More](#)

Parameter	Value
Product Name	YIPF7 Rabbit pAb
Remarks/Alias	FinGER9; YIPF7
Species	Human
GeneID (Human)	285525
GeneID	285525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human YIPF7 (NP_872398.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q8N8F6
Protein Weight	31kDa
Shipping	Ice bag