

Anti-Ykl-40/Chi3L1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P36222

Item Number: CM0003455



Introduction

YKL-40/CHI3L1 Rabbit mAb (recombinant) for WB, IF/ICC, IF-P, ELISA. Specific to human YKL-40, cross-reacts with mouse, rat. Ideal for inflammation, tissue remodeling, and biomarker research.

[Learn More](#)

Parameter	Value
Product Name	YKL-40/CHI3L1 Rabbit mAb
Remarks/Alias	GP39; ASRT7; GP-39; YK-40; YKL40; CGP-39; YKL-40; YYL-40; HC-gp39; HCGP-3P; hCGP-39; YKL-40/CHI3L1
Species	Human
GeneID	1116 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 284-383 of human YKL-40/CHI3L1 (NP_001267.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36222
Protein Weight	43 kDa
Shipping	Ice bag