


Anti-Ykl-40/Chi3L1 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q61362

Item Number: CM0011075



Introduction

Anti-YKL-40/CHI3L1 rabbit monoclonal antibody validated for Mouse samples in WB, IP, IHC-P, ELISA. Targets chitinase-3-like protein 1 (CHI3L1/BRP39), involved in tissue remodeling and Th2 inflammation.

[Learn More](#)

Parameter	Value
Product Name	YKL-40/CHI3L1 Rabbit mAb
Remarks/Alias	Gp39; Brp39; Chi1
Species	Mouse
GeneID (Human)	1116
GeneID (Mouse)	12654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-389 of Mouse Chi3l1 (NP_031721.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q61362
Protein Weight	43 kDa/44 kDa
Shipping	Ice bag