

Anti-Ykl-40/Chi3L1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36222

Item Number: CM0030030



Introduction

YKL-40/CHI3L1 Rabbit pAb validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat CHI3L1, a 43kDa lectin involved in inflammation and tissue remodeling. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	YKL-40/CHI3L1 Rabbit pAb
Remarks/Alias	GP39; ASRT7; GP-39; YK-40; YKL40; CGP-39; YKL-40; YYL-40; HC-gp39; HCGP-3P; hCGP-39; YKL-40/CHI3L1
Species	Human
GeneID (Human)	1116
GeneID	1116
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-383 of YKL-40/CHI3L1 (NP_001267.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36222
Protein Weight	43kDa
Shipping	Ice bag