

Anti-Yme1L1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ta2

Item Number: CM0016611

Introduction

YME1L1 Rabbit Polyclonal Antibody designed for WB and ELISA, recognizing human, mouse, and rat YME1L1. Ideal for investigating mitochondrial proteostasis and morphology. This antibody targets the N-terminal region of human YME1L1.

[Learn More](#)

Parameter	Value
Product Name	YME1L1 Rabbit pAb
Remarks/Alias	FTSH; MEG4; PAMP; OPA11; YME1L; YME1L1
Species	Human
GeneID (Human)	10730
GeneID	10730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human YME1L1 (NP_055078.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96TA2
Protein Weight	86kDa
Shipping	Ice bag