


Anti-Yod1 Rabbit Monoclonal Antibody (Wb, Elisa) - Q5Vvq6

Item Number: CM0006289



Introduction

Anti-YOD1 Rabbit Monoclonal Antibody validated for Western blotting (WB) and ELISA. Detects human YOD1, a deubiquitinase involved in ERAD and autophagy. Rabbit-derived monoclonal antibody.

[Learn More](#)

Parameter	Value
Product Name	YOD1 Rabbit mAb
Remarks/Alias	DUBA8; OTUD2; PRO0907; YOD1
Species	Human
GeneID	55432
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 249-348 of human YOD1 (NP_061036.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5VVQ6
Protein Weight	38 kDa
Shipping	Ice bag