

Anti-Ypel3 Polyclonal Antibody For Wb, Elisa - P61236

Item Number: CM0024345

Introduction

Rabbit polyclonal antibody against human YPEL3, validated in Western blot and ELISA, cross-reactive with mouse and rat. Ideal for apoptosis and cell proliferation studies.

[Learn More](#)

Parameter	Specification
Product Name	YPEL3 Rabbit pAb
Remarks/Alias	YPEL3
Species	Human
Gene ID (Human)	83719
Gene ID	83719
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human YPEL3 (NP_001138996.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P61236
Protein Weight	14kDa
Shipping	Ice bag