


Anti-Ypel5 Polyclonal Antibody For Wb, Elisa - P62699

Item Number: CM0013651



Introduction

YPEL5 Rabbit Polyclonal Antibody validated for WB and ELISA, reactive with Human, Mouse, Rat. Targets yippee-like 5, a CTLH E3 ligase component involved in cell proliferation. Molecular weight: 14kDa.

[Learn More](#)

Parameter	Value
Product Name	YPEL5 Rabbit pAb
Remarks/Alias	CGI-127; YPEL5
Species	Human
GeneID (Human)	51646
GeneID	51646
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62699
Protein Weight	14kDa
Shipping	Ice bag