

Anti-Ythdc1 Monoclonal Antibody (Ko Validated) For Wb - E9Q5K9

Item Number: CM0004070

Introduction

KO Validated rabbit monoclonal antibody against mouse Ythdc1 (E9Q5K9), for Western blot and ELISA applications. Recombinant protein immunogen, 85 kDa. Cross-reacts with human YTHDC1.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Ythdc1 Rabbit mAb
Remarks/Alias	mKIAA1966; A730098D12Rik; Mouse Ythdc1
Species	Mouse
Gene ID (Human)	91746
Gene ID	231386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 597-736 of mouse Ythdc1 (NP_808348.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	E9Q5K9
Protein Weight	85kDa
Shipping	Ice bag