


Anti-Ythdc1 Polyclonal Antibody For Wb And Elisa - Q96Mu7

Item Number: CM0029631



Introduction

Rabbit polyclonal antibody against human YTHDC1, validated for WB and ELISA, cross-reactive with mouse and rat. Targets the N6-methyladenosine reader protein involved in alternative splicing and RNA metabolism.

[Learn More](#)

Parameter	Value
Product Name	YTHDC1 Rabbit pAb
Remarks/Alias	YT521; YT521-B; YTHDC1
Species	Human
Gene ID (Human)	91746
Gene ID	91746
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96MU7
Protein Weight	85kDa
Shipping	Ice bag