

Anti-Ythdf1 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Byj9

Item Number: CM0015978



Introduction

This rabbit anti-YTHDF1 polyclonal antibody detects human, mouse, and rat YTHDF1. Validated for Western blot, IF/ICC, IP, and ELISA. Useful for m6A RNA modification, mRNA metabolism, and epitranscriptomics research.

[Learn More](#)

Parameter	Value
Product Name	YTHDF1 Rabbit pAb
Remarks/Alias	DF1; C20orf21; YTHDF1
Species	Human
GeneID (Human)	54915
GeneID	54915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human YTHDF1 (NP_060268.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYJ9
Protein Weight	61kDa
Shipping	Ice bag