

Anti-Ythdf1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Byj9

Item Number: CM0021903



Introduction

Rabbit polyclonal antibody against YTHDF1, recognizing human, mouse, and rat. Validated for Western blot, IHC-P, and ELISA, with a predicted molecular weight of 61kDa. Ideal for m6A RNA modification and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	YTHDF1 Rabbit pAb
Remarks/Alias	DF1; C20orf21; YTHDF1
Species	Human
GeneID (Human)	54915
GeneID	54915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-390 of human YTHDF1
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYJ9
Protein Weight	61kDa
Shipping	Ice bag