

Anti-Ythdf1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Byj9

Item Number: CM0006467



Introduction

YTHDF1 Rabbit mAb for WB, IP, ELISA – specifically recognizes m6A-modified mRNAs, regulating stability and translation. Validated in human, mouse, rat. Ideal for epitranscriptomics, cancer, and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	YTHDF1 Rabbit mAb
Remarks/Alias	YTHDF1; C20orf21; YTH domain-containing family protein 1
Species	Human
GeneID (Human)	54915
GeneID	54915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240aa of human YTHDF1 (NP_060268.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYJ9
Protein Weight	20kDa/60kDa
Shipping	Ice bag