

Anti-Ythdf2 (Ko Validated) Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - Q9Y5A9

Item Number: CM0010065



Introduction

KO validated rabbit monoclonal antibody against YTHDF2, an m6A reader protein. Recognizes human, mouse, rat. Suitable for WB, IP, IF/ICC, ELISA. Ideal for studying RNA metabolism, epitranscriptomics, and viral infection.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] YTHDF2 Rabbit mAb
Remarks/Alias	DF2; CAHL; HGRG8; NY-REN-2
Species	Human
GeneID (Human)	51441
GeneID	51441
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 339-565 of human YTHDF2 (NP_057342.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5A9
Protein Weight	62 kDa/57 kDa
Shipping	Ice bag