

Anti-Ythdf2 Rabbit Polyclonal Antibody For Wb, Ihc, If - Q9Y5A9

Item Number: CM0021246



Introduction

Rabbit polyclonal antibody against YTHDF2, an m6A reader protein. Validated for WB, IHC, IF, and ELISA. Reacts with human, mouse, and rat. Ideal for epitranscriptomics and RNA biology research.

[Learn More](#)

Parameter	Value
Product Name	YTHDF2 Rabbit pAb
Remarks/Alias	DF2; CAHL; HGRG8; NY-REN-2; YTHDF2
Species	Human
GeneID (human)	51441
GeneID	51441
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human YTHDF2 (NP_057342.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5A9
Protein Weight	62kDa
Shipping	Ice bag