

Anti-Ythdf3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q7Z739

Item Number: CM0027486



Introduction

Rabbit polyclonal antibody against YTHDF3, an m6A reader regulating mRNA stability. Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. Ideal for epitranscriptomics research.

[Learn More](#)

Parameter	Value
Product Name	YTHDF3 Rabbit pAb
Remarks/Alias	DF3; YTHDF3
Species	Human
GeneID (Human)	253943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human YTHDF3 (NP_689971.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z739
Protein Weight	64kDa
Shipping	Ice bag