

Anti-Yucca1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Szy8

Item Number: CM0018921

Introduction

YUCCA1 Rabbit Polyclonal Antibody detects Arabidopsis thaliana YUCCA1 protein (~46 kDa) by WB and ELISA. YUCCA1 is a flavin-containing monooxygenase involved in auxin biosynthesis and plant development. Ideal for plant research.

[Learn More](#)

Parameter	Value
Product Name	YUCCA1 Rabbit pAb
Remarks/Alias	L23H3.20; L23H3_20; YUC; YUCCA; YUCCA 1; YUCCA1
Species	Arabidopsis thaliana
GeneID	829389
Immunogen	Recombinant protein corresponding to amino acids 71-221 of Arabidopsis thaliana YUCCA1 (NP_194980.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q9SZY8
Protein Weight	46kDa
Shipping	Ice bag