


Anti-Ywhah Polyclonal Antibody For Wb, If/Icc, Elisa - Q04917

Item Number: CM0030208



Introduction

YWHAH Rabbit pAb detects 14-3-3 eta (28 kDa) via WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for brain signaling studies and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	YWHAH Rabbit pAb
Remarks/Alias	YWHA1; YWHAH
Species	Human
Gene ID (Human)	7533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human YWHAH (NP_003396.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04917
Protein Weight	28kDa
Shipping	Ice bag