


Anti-Ywhaz (14-3-3 Zeta/Delta) Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P63104

Item Number: CM0030341



Introduction

High-quality rabbit polyclonal antibody against human YWHAZ (14-3-3 zeta/delta). Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for signaling pathway research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	YWHAZ Rabbit pAb
Remarks/Alias	HEL4; YWHAD; KCIP-1; HEL-S-3; POPCHAS; HEL-S-93; 14-3-3-zeta; YWHAZ
Species	Human
GeneID (Human)	7534
GeneID	7534
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 146-245 of human YWHAZ (NP_003397.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63104
Protein Weight	28kDa
Shipping	Ice bag