

Anti-Yy1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa, Chip - P25490

Item Number: CM0002403



Introduction

This rabbit monoclonal antibody against human YY1 (P25490) is validated for WB, IF/ICC, ELISA, and ChIP. Cross-reacts with mouse and rat. Ideal for studying transcriptional regulation and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	YY1 Rabbit mAb
Remarks/Alias	DELTA; NF-E1; UCRBP; GADEV5; INO80S; YIN-YANG-1; YY1
Species	Human
GeneID (Human)	7528
GeneID	7528
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 200-300 of human YY1 (P25490).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA, CHIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P25490
Protein Weight	45kDa
Shipping	Ice bag