

Anti-Yy1Ap1 Polyclonal Antibody For Wb And Elisa - Q9H869

Item Number: CM0015053

Introduction

Rabbit polyclonal antibody targeting human YY1AP1 (HCCA2). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for chromatin remodeling and liver cancer research.

[Learn More](#)

Parameter	Value
Product Name	YY1AP1 Rabbit pAb
Remarks/Alias	GRNG; HCCA1; HCCA2; YY1AP; YY1AP1
Species	Human
Gene ID (Human)	55249
Gene ID	55249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human YY1AP1 (NP_060723.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H869
Protein Weight	88kDa
Shipping	Ice bag