


Anti-Zap70 Polyclonal Antibody For Wb, If/Icc, Elisa - P43403

Item Number: CM0019750



Introduction

Rabbit polyclonal antibody against human ZAP70 (immunogen: aa 1-300). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for T-cell signaling and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	ZAP70 Rabbit pAb
Remarks/Alias	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70; ZAP70
Species	Human
GeneID (human)	7535
GeneID	7535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ZAP70 (NP_001070.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43403
Protein Weight	70kDa
Shipping	Ice bag