


Anti-Zap70 Rabbit Monoclonal Antibody For Wb, Elisa - P43403

Item Number: CM0009585



Introduction

Recombinant rabbit monoclonal anti-ZAP70 antibody (70kDa). Validated for WB and ELISA; cross-reacts with mouse. Ideal for T-cell signaling, lymphocyte development, and adaptive immunity research.

[Learn More](#)

Parameter	Specification
Product Name	ZAP70 Rabbit mAb
Remarks/Alias	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70; ZAP70
Species	Human
GeneID (Human)	7535
GeneID	7535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-619 of human ZAP70 (NP_001070.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P43403
Protein Weight	70kDa
Shipping	Ice bag