


Anti-Zbed1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O96006

Item Number: CM0029032



Introduction

Rabbit polyclonal antibody against human ZBED1, validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for transcription and cell proliferation research; detects endogenous ZBED1 at 78 kDa.

[Learn More](#)

Parameter	Value
Product Name	ZBED1 Rabbit pAb
Remarks/Alias	ALTE; DREF; TRAMP; hDREF; ZBED1
Species	Human
GeneID (human)	9189
GeneID	9189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZBED1 (NP_001164607.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O96006
Protein Weight	78kDa
Shipping	Ice bag