

Anti-Zbp1 Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q9Qy24

Item Number: CM0026603



Introduction

Rabbit polyclonal antibody against mouse ZBP1 (DAI/DLM-1) for WB, IF/ICC, IF-P, IHC-P, and ELISA. Recognizes Z-DNA-binding protein 1, a key innate immune sensor triggering PANoptosis.

[Learn More](#)

Parameter	Value
Product Name	ZBP1 Rabbit pAb
Remarks/Alias	Dai; Dlm1; mZaDLM; 2010010H03Rik
Species	Mouse
GeneID (Human)	81030
GeneID	58203
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of mouse ZBP1 (NP_067369.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q9QY24
Protein Weight	44 kDa/21 kDa
Shipping	Ice bag