

Anti-Zbp1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q9H171

Item Number: CM0007739



Introduction

Rabbit monoclonal antibody targeting human ZBP1 (DAI/DLM-1), validated for WB, IP, IF/ICC, IHC-P, ELISA. Detects 46kDa Z-DNA/RNA binding protein, a key innate immune sensor of viral Z-RNA triggering PANoptosis.

[Learn More](#)

Parameter	Value
Product Name	ZBP1 Rabbit mAb
Remarks/Alias	DAI; DLM1; DLM-1; C20orf183
Species	Human
GeneID (Human)	81030
GeneID	81030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human ZBP1 (NP_110403.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9H171
Protein Weight	46kDa
Shipping	Ice bag