

Anti-Zbp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H171

Item Number: CM0014782



Introduction

High-quality rabbit polyclonal antibody against ZBP1 (DAI/DLM-1), validated for Western blot and ELISA, reactive with human, mouse, and rat samples. Ideal for studying innate immunity and PANoptosis pathways. - Q9H171

[Learn More](#)

Parameter	Value
Product Name	ZBP1 Rabbit pAb
Remarks/Alias	DAI; DLM1; DLM-1; C20orf183; ZBP1
Species	Human
Gene ID (Human)	81030
Gene ID	81030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human ZBP1 (NP_110403.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H171
Protein Weight	46kDa
Shipping	Ice bag