


Anti-Zbtb17 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13105

Item Number: CM0026849



Introduction

Rabbit polyclonal antibody against ZBTB17 (MIZ-1), a transcription factor critical for lymphocyte development and cell cycle regulation. Validated for WB, IHC-P, ELISA, and reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	ZBTB17 Rabbit pAb
Remarks/Alias	MIZ-1; ZNF60; ZNF151; pHZ-67; ZBTB17
Species	Human
GeneID (Human)	7709
GeneID	7709
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-360 of human ZBTB17 (NP_003434.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13105
Protein Weight	88kDa
Shipping	Ice bag