

Anti-Zbtb6 Rabbit Polyclonal Antibody For Wb, Elisa - Q15916

Item Number: CM0015424

Introduction

Rabbit polyclonal antibody against ZBTB6 (ZID, ZNF482) for WB and ELISA. Detects human, mouse, and rat ZBTB6, a zinc finger transcription factor broadly expressed with highest levels in brain.

[Learn More](#)

Parameter	Value
Product Name	ZBTB6 Rabbit pAb
Remarks/Alias	ZID; ZNF482; ZBTB6
Species	Human
GeneID (human)	10773
GeneID	10773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ZBTB6 (NP_006617.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15916
Protein Weight	48kDa
Shipping	Ice bag